

ADA Dropln Sinks

Innovations for Accessibility



SECTION

DROP-IN SINK

SINK	WEIGHT
ADA26	39 LBS

Phone: 512/352-5591
FAX: 512/352-9523

ADA DropIn® Sinks




ADA DropIn® Sink

-sink supports. Top mounting also elevates the already shallow sink bottom by about 1". All outside edges and corners have been rounded for the user's comfort and protection and all ADA DropIn sinks feature space saving rear corner drains to move plumbing to the back of the cabinet.


LabTops' Tiered ADA DropIn Sink features a graduated depth basin to allow wheelchair access to sinks with full depth basins. The 5" deep front tier can also be used to place small items or apparatus for drying while rinsing or cleaning large items in the 11" deep end.

DropIn Sink Systems can be provided by Laboratory Tops to allow installation of regular depth sinks that can be changed out to deeper basin sinks or to a shallow basin ADA DropIn sink as needed. Two DropIn Sink Systems include ADA DropIns* (please refer to our Product Bulletin 907 for further information on these systems).



Tiered ADA DropIn® Sink

For further ADA information please reference the U.S. Department of Justice's Americans with Disabilities Act home page <http://www.usdoj.gov/crt/ada/adahom1.htm> (Note: Some states have additional requirements that supersede federal guidelines).



ADA DropIn® Sink Sizes

Sink Styles	Inside Dimensions			Outside Dimensions			Inside Dimensions	
	A	B	C	D	E	F	G	H
A05	14.0	10.0	5.0	15.6	11.6	5.8	3.5	3.5
A07	14.0	14.0	5.0	15.6	15.6	5.8	3.5	3.5
A25*	18.0	15.0	5.0	19.6	16.6	5.8	3.5	3.5
A26*	18.0	15.0	5/11	19.6	16.6	12	5.5	4.5
A55*	25.0	15.0	4.8	26.6	16.6	5.8	3.5	3.5

ADA DropIn Sinks are available through laboratory casework suppliers or directly from Laboratory Tops, Inc. If you need further information about these or any other quality LabTops products, please contact our customer service representatives at:

